

Core Board Teaching Guide

A visual map of the exact 12 x 7 Full View board, its color system, and repeatable paths for teaching language.

No account	No ads	No cloud profile	Works offline after loading
------------	--------	------------------	-----------------------------

I	we	want	like	go	help	is	on	in	up	feel	here
you	they	play	eat	drink	look	come	give	do	can	talk	again
he	she	make	get	put	open	turn	read	sit	walk	run	down
it	sleep	need	that	love	happy	not	out	off	sad	and	the
yes	more	good	big	hot	some	all done	little	bad	no	stop	cold
my	please	thank you	sorry	what	where	who	when	why	how	to	grammar
CLEAR	food	people	feelings	places	animals	things	colors	body	clothes	my words	ABC

Motor planning rule: a Level 1 word stays in the same physical location when its related word page opens. Teach the movement, not just the picture.

Why the board is organized this way

- Rows 1-4 hold pronouns, actions, grammar, and descriptions used across many situations.
- Row 5 adds responses and describing words.
- Row 6 holds social words, question words, to, and the grammar doorway.
- Row 7 contains category doorways. The category is a path to vocabulary, not a spoken sentence.

Color key

	Pronouns	I, you, he, she, we, they
	Actions	want, go, make, help, feel
	Descriptions	happy, big, hot, little
	People & things	school, food, people, objects
	Social words	yes, no, please, stop, sorry
	Grammar	to, the, and, in, where, why
	Folders	food, people, places, body, ABC

Teach color as a clue, never as the answer

Say “actions are green, so look in the green area” while still naming and modeling the word. A communicator may use location, color, picture, print, or all four.

Useful Level 1 paths

Goal	Model	Why it matters
Request	I want more	Pronoun + action + description
Go somewhere	I want to go to school	Includes the grammar word to
Reject	I do not want that	A complete refusal
Ask	where is my phone	Question + grammar + possession
Comment	I like that	Communication beyond requesting
Health	I feel sick and need help	Describes and requests support
Repair	not that, I want this	Corrects a misunderstanding

Modeling routes

Do not move words to make them easier. Repeat the same route until it becomes familiar.

Message	Route to model
go to school	Home: go -> Home: to -> Home: school / Places: school
I need bathroom	Home: I -> Home: need -> Places/Body vocabulary: bathroom
open the door	Home: open -> Home: the -> Things: door
my tummy hurts	Home: my -> Body: tummy -> Health/Body: hurt
what happened	Home: what -> action page: happened
I want a different game	Home: I -> want -> Grammar: a -> Descriptions: different -> Things: game

Use Find as a teaching tool

- Type a word and choose a result. The app highlights each press from Home in order.
- Watch the path, then return Home and model the same path yourself.
- If no exact word appears, add a useful personal word to the most meaningful folder.

Simple View

Simple View uses a 6 x 8 layout with larger targets for phones and smaller screens. It is a complete communication view, and repeated words keep stable Simple View positions.

Communication is successful when the message is understood - even if the sentence is incomplete, unconventional, or made with several methods.